

Appendix I Specifications

A1.1 Recommended Hardware Configuration

The minimum requirements of hardware configuration for the CMS-LITE are shown as below.

Components	Requirements											
System	Meet the IEC/EN control requirements for ITE device											
PC workstation	CPU: Dual Intel Core i3 3.0 Ghz or above Memory: 4 G or above Hard disk: 500 GB or above Display interface: 2 LAN port: 1 or above USB port: more than one OS: Windows 10 (Above the version 1507.10240) Pro 32/64 bit											
Keyboard	PS/2 or USB keyboard with CE mark											
Mouse	PS/2 or USB mouse with CE mark											
Display	Specifications: <table border="1" style="margin-left: 20px;"> <thead> <tr> <th>Dimensions (inch)</th> <th>Resolution (pixel)</th> </tr> </thead> <tbody> <tr> <td>19 (widescreen)</td> <td>1440X900</td> </tr> <tr> <td>*19 (regular-screen)</td> <td>1280X1024</td> </tr> <tr> <td>22 (widescreen)</td> <td>1920X1080</td> </tr> <tr> <td colspan="2"><i>*Recommended</i></td> </tr> </tbody> </table>		Dimensions (inch)	Resolution (pixel)	19 (widescreen)	1440X900	*19 (regular-screen)	1280X1024	22 (widescreen)	1920X1080	<i>*Recommended</i>	
Dimensions (inch)	Resolution (pixel)											
19 (widescreen)	1440X900											
*19 (regular-screen)	1280X1024											
22 (widescreen)	1920X1080											
<i>*Recommended</i>												
Printer	LaserJet											
UPS	1000 W											
Network device specifications	Structure	Ethernet 802.3										
	Device	Network switch										
	Transmission rate	10 M, 100 M										
	Transmitted information	Waveforms, parameters and alarms of all networked bedside monitors										
	Compatible monitors	Patient monitors complying with EDAN network protocol										