

Item	Client Specification	Server Specification
System	Satisfies the IEC/EN control requirements	Satisfies the IEC/EN control requirements
PC workstation	Microcomputer Model: POC-W212C	Microcomputer Model: POC-W212C
Keyboard	PS/2 or USB keyboard with CE-mark	PS/2 or USB keyboard with CE-mark
Mouse	PS/2 mouse with CE-mark	PS/2 mouse with CE-mark
Display	Resolution: 1366×768 or above	Resolution: 1024×768 or above
Printer	HP LaserJet series with CE Mark Parallel port or USB port Printer performance specifications: Resolution: 600DPI * 600DPI Memory: 1 MB Interface: IEEE 1284 Parallel port or USB port Speed: 6 pages / min. Paper size: A4/Letter	/
UPS	/	1000 W with CE Mark (required when the electricity fluctuates abnormally)
Speaker	The speaker of POC-W212C	/

**CAUTION**

- 1 Make sure the computer hardware meets the installing and running requirements of the software, otherwise the software may not work properly.
- 2 If the host crashes, you need to restart the computer.
- 3 Install anti-virus software on your computer and run it regularly in case the system is subject to viral attacks.
- 4 The hardware specifications require the use of PC that complies with IEC requirements for ITE equipment.
- 5 The configuration mentioned above is not permanent. THE MANUFACTURER preserves the right to change and update system settings.
- 6 Need to meet the recommended resolution or higher, in order to achieve monitoring data to see the effect.

**1.3 Operating System Requirements**

The OS of the server computer can be:

- ◆ Microsoft Windows 7 SP1 (32/64 bit)
- ◆ Microsoft Windows 8(32/64 bit)
- ◆ Microsoft Windows 8.1(32/64 bit)
- ◆ Microsoft Windows 10(32/64 bit), 1809 version
- ◆ Windows Server 2008
- ◆ Windows Server 2012

The OS of the client computer can be:

- ◆ Microsoft Windows 7 SP1 (32/64 bit)
- ◆ Microsoft Windows 8(32/64 bit)
- ◆ Microsoft Windows 8.1(32/64 bit)
- ◆ Microsoft Windows 10(32/64 bit), 1809 version